

Enhanced Streptavidin-HRP

LABELS & CONJUGATES

Streptavidin, poly-HRP linker Cat. no. 4740 (1 mg HRP/mL)

Product Characteristics

Enhanced Streptavidin-HRP is a HRP-conjugate specific for biotinylated antibodies and is characterized by an approximate enhanced molar ratio of >4 HRP enzymes covalently coupled per streptavidin molecule. Enhanced Streptavidin-HRP conjugate enables at least a 4-fold increase in sensitivity in ELISA compared to conventional conjugate and offers pg/fg sensitivity levels in immunoassays, when only minimal amounts of specific antibodies or antigens are available:

- Signal amplification
- Ultrasensitive detection
- Faster assay time
- Recommended when only minimal amounts of specific antibody or antigen are available
- Can be used with any biotin labeled probe

Composition & Properties

Specificity	: Biotinylated probes
Content	: Stabilized in 50% TRIS Buffer and 50% HRP Conjugate Stabilizer
Activity	: 90 - 110%
Recommended Dilution	: 1:15.000 to 1:30.000

Tips & Tricks

- Enhanced Streptavidin-HRP can successfully be diluted and stored in Antibody Enhancer (cat. no. 4660). (Not to be used for dilution of Ab's directed against rabbit immunoglobulins).
- Avoid batch-to-batch variations and day-to-day variation by preparing a stock solution in your own dilution.
- Fine-tuning of the polymer backbone is necessary for Immunoblotting, IHC & DNA probe detection.

Handling & Storage

Store solution at 2-8 °C.



Our Immunoassay Solutions are eco-friendly, creating a healthy work environment and preserving natural resources, while helping our customers fulfill significant regulatory requirements. We call it ECO-TEK.

Gentaur Molecular Products BVBA Address: Voortstraat 49, 1910 Kampenhout, Belgium T: 0032 16 58 90 45 | E: info@gentaur.com Websites: www.gentaur.com | www.maxanim.com