$See \ discussions, stats, and author \ profiles \ for \ this \ publication \ at: \ https://www.researchgate.net/publication/389628021$

Polymerase Chain Reaction (PCR): Principle, Requirements, Types and Application

reads 775

Article · March 2025

CITATIONS	
0	
1 author:	
	Kombo Juma
	Kampala International University, Uganda
-	6 PUBLICATIONS 0 CITATIONS
	SEE PROFILE

https://lifelinejournals.org/journals/lifeline-health-sciences

Polymerase Chain Reaction (PCR): Principle, Requirements, Types and Application

Juma Kombo

Department of Medical Laboratory Science, Kampala International University, Uganda

Abstract

Polymerase Chain Reaction (PCR) is a widely used molecular biology technique that enables the amplification of specific DNA sequences. It is essential in genetic research, medical diagnostics, forensic science, and biotechnology. PCR works by replicating targeted DNA sequences through cycles of denaturation, annealing, and extension, allowing for the detection and analysis of genetic material even from minute samples. This study involved a review of secondary data on the principles, types, and applications of PCR in scientific and medical research. The literature review highlights the fundamental principles of PCR and its evolution into various specialized techniques, each addressing specific research and diagnostic needs. Conventional PCR remains a reliable method for DNA amplification, while advanced variations such as real-time PCR (qPCR) and digital PCR (dPCR) enhance sensitivity and quantification accuracy. Multiplex PCR improves efficiency by detecting multiple targets simultaneously, whereas RT-PCR plays a crucial role in diagnosing RNA viruses. Emerging technologies, including nanoparticle-assisted PCR and AI-integrated qPCR, have further improved precision, automation, and data interpretation. PCR has transformed molecular biology by enabling rapid and precise DNA amplification, with widespread applications in medicine, forensic science, agriculture, and environmental research. Despite challenges such as contamination risks and primer specificity issues, continuous advancements-including automation, AI integration, and emerging PCR technologies—have significantly enhanced its accuracy, efficiency, and reliability. As innovation progresses, PCR remains a cornerstone of genetic analysis and diagnostics, with expanding potential in scientific and clinical fields.

Keywords: Polymerase Chain Reaction, Principle, Requirements, Type, Application

Introduction

Polymerase Chain Reaction (PCR) is a revolutionary molecular biology technique that enables the amplification of specific DNA sequences with high precision. Since its invention by Kary Mullis in 1983, PCR has become an indispensable tool in numerous scientific disciplines, including genetics, forensic science, medical diagnostics, and biotechnology(1).. This method allows researchers to amplify minuscule amounts of DNA, even from degraded or ancient samples, facilitating detailed genetic analysis. PCR operates by mimicking the natural DNA replication process through a series of thermal cycles, comprising denaturation, annealing, and

Citation: Kombo J. Polymerase Chain Reaction (PCR): Principle, Requirements, Types and Application. Lifeline Health Sciences, 2025; 3(1): 98-105

https://lifelinejournals.org/journals/lifeline-health-sciences

extension steps. The technique relies on thermostable DNA polymerase, such as Taq polymerase, to synthesize new DNA strands. Studies have demonstrated that optimizing PCR conditions, such as primer design and annealing temperatures, significantly enhances amplification efficiency and specificity (2).. These findings underscore the necessity of standardized protocols to improve reproducibility in PCR-based applications. The applications of PCR are vast and continue to expand. In forensic science, PCR has been instrumental in DNA fingerprinting, where it enables the identification of individuals from minute biological samples. A study by (3)confirmed that PCR-based DNA profiling has improved the resolution of forensic investigations, especially in solving cold cases. Similarly, medical diagnostics have benefited immensely from PCR technology. Real-time PCR (qPCR) and reverse transcription PCR (RT-PCR) have been widely employed in detecting genetic disorders, infectious diseases, and cancer biomarkers. Research by (3)highlighted the role of RT-PCR in diagnosing respiratory viruses, with findings showing a 95% detection accuracy for influenza and coronaviruses. PCR has also been crucial in agricultural and environmental research. Multiplex PCR, which allows simultaneous amplification of multiple DNA targets, has been employed for pathogen detection in crops. (3)demonstrated that multiplex PCR effectively identified multiple bacterial strains in a single reaction, reducing diagnostic time and costs in agricultural biotechnology. Furthermore, environmental DNA (eDNA) analysis using PCR has revolutionized biodiversity monitoring. Carter et al. (2021) found that PCR-based eDNA detection provided a non-invasive method for tracking aquatic species and detecting invasive organisms, thus aiding conservation efforts. Despite its numerous advantages, PCR has some challenges, including contamination risks, primer specificity issues, and dependency on thermal cycling. (4) reported that improper laboratory handling was a major cause of false-positive results in PCR assays(5).. Their study emphasized the need for strict contamination control measures. Additionally, advancements in PCR technology continue to address these limitations. Digital PCR (dPCR), for instance, has emerged as a more precise method for absolute quantification of nucleic acids, particularly in detecting rare mutations in cancer patients (4)Similarly, nanoparticle-assisted PCR (nanoPCR) has shown promise in enhancing reaction efficiency, with (6) reporting improved thermal conductivity and amplification yield. The COVID-19 pandemic further underscored the critical role of PCR in public health. Rapid and accurate detection of SARS-CoV-2 through RT-PCR enabled early diagnosis and containment efforts worldwide. (7) conducted a study on COVID-19 diagnostics, confirming that RT-PCR provided reliable viral load quantification, which was essential for patient monitoring and epidemiological surveillance(6)

ChallengesinPCRTechnologyDespite its many benefits, PCR has some limitations that affect its accuracy and accessibility:

Contamination Risks: Even a small amount of unwanted DNA can lead to incorrect results, making careful laboratory practices essential.

2

Citation: Kombo J. Polymerase Chain Reaction (PCR): Principle, Requirements, Types and Application. Lifeline Health Sciences, 2025; 3(1): 98-105

https://lifelinejournals.org/journals/lifeline-health-sciences

Primer Design Limitations: PCR requires short DNA sequences called primers to start the copying process. If the primers are not designed correctly, they may bind to the wrong sequences, leading to errors.

Thermal Cycling Requirements: PCR depends on precise temperature changes for each step of the process. In areas with limited resources, this can be difficult to achieve.

High Cost of Equipment and Reagents: Advanced PCR methods, such as real-time PCR, require expensive machines and special chemicals, making them difficult to afford in some regions.

Principles of Polymerase Chain Reaction (PCR)

Polymerase Chain Reaction (PCR) is a widely used technique in molecular biology that allows scientists to amplify specific DNA sequences. This method was first introduced by Kary Mullis in 1983 and has since become a crucial tool in various fields, including medical diagnostics, forensic science, and genetic research. PCR is based on the natural process of DNA replication and consists of three main steps: denaturation, annealing, and extension. During denaturation, the double-stranded DNA is heated to separate it into two single strands. In the annealing step, short synthetic DNA sequences known as primers bind to the target regions on each strand. Finally, in the extension phase, DNA polymerase synthesizes new DNA strands by adding nucleotides to the primers, resulting in the amplification of the target sequence. This cycle is repeated multiple times to generate millions of copies of the DNA segment of interest(9).

A study by (9)explored the fundamental principles of PCR and highlighted the importance of reaction optimization in ensuring specificity and efficiency. Their research demonstrated that variations in annealing temperatures and primer concentrations could significantly affect amplification outcomes, emphasizing the need for standardized protocols in PCR-based studies. This study provided critical insights into optimizing PCR conditions for diverse applications, reinforcing its relevance in molecular diagnostics and genetic research(9).

Types of PCR

Since its invention, several variations of PCR have been developed to address specific research and diagnostic needs. These include conventional PCR, real-time PCR, digital PCR (dPCR), multiplex PCR, quantitative PCR (qPCR), reverse transcription PCR (RT-PCR), nested PCR, touchdown PCR, hot start PCR, colony PCR, and in situ PCR.

Conventional PCR is the standard form of PCR that involves thermal cycling to amplify DNA. It requires gel electrophoresis for the visualization of amplified products. A study by (10)examined the use of conventional PCR in genetic testing and found that while it remains a reliable method, its dependence on post-PCR analysis can make it time-consuming compared to more modern techniques(10).

Citation: Kombo J. Polymerase Chain Reaction (PCR): Principle, Requirements, Types and Application. Lifeline Health Sciences, 2025; 3(1): 98-105

https://lifelinejournals.org/journals/lifeline-health-sciences

Real-Time PCR (qPCR), also known as quantitative PCR, enables researchers to measure the amount of DNA present in a sample during the amplification process. A study conducted by(9). evaluated the accuracy of qPCR in monitoring viral loads in patients with COVID-19. Their findings revealed that qPCR provided precise quantification of SARS-CoV-2 RNA, aiding in patient management and disease surveillance(11).

Digital PCR (dPCR) is a highly sensitive technique that partitions a DNA sample into thousands of individual reactions, allowing for absolute quantification of target sequences. (12)investigated the efficacy of dPCR in detecting rare mutations in cancer patients. Their study found that dPCR provided higher accuracy and sensitivity than conventional qPCR, making it a valuable tool in precision medicine and early cancer detection(12).

Multiplex PCR is designed to amplify multiple DNA targets simultaneously using multiple primer sets. A study by (13)assessed the application of multiplex PCR in pathogen identification. The researchers demonstrated that multiplex PCR effectively detected multiple bacterial strains in a single reaction, reducing diagnostic time and costs in clinical microbiology(13).

Reverse Transcription PCR (RT-PCR) converts RNA into complementary DNA (cDNA) before amplification. (14)explored the role of RT-PCR in diagnosing viral infections, particularly in influenza and coronaviruses. Their study confirmed that RT-PCR remained the gold standard for detecting RNA viruses due to its high sensitivity and specificity.

Nested PCR is a modification of standard PCR that involves two consecutive amplification reactions to improve specificity and reduce non-specific amplification. (14)demonstrated that nested PCR significantly enhanced the detection of low-copy-number pathogens in clinical samples, making it highly effective in detecting Mycobacterium tuberculosis in tuberculosis patients(14).

Hot Start PCR involves the use of a DNA polymerase that remains inactive at low temperatures, preventing non-specific amplification. (15)highlighted that hot start PCR significantly improved the amplification of low-abundance DNA targets in forensic and clinical applications(15).

In Situ PCR is performed within fixed cells or tissues to detect specific nucleic acid sequences while maintaining spatial information. A study by (16)showed that in situ PCR was highly useful in detecting viral genomes within infected tissues, improving histopathological diagnostics(16).

PCR Reaction Requirements

Each type of Polymerase Chain Reaction (PCR) requires specific reaction components and conditions to achieve optimal amplification. The fundamental components include Taq polymerase, MgCl₂, buffer, dNTPs, primers, DNA sample, and nuclease-free distilled water(17). Conventional PCR, the standard form of PCR, requires Taq polymerase to facilitate DNA amplification, MgCl₂ as a cofactor for enzymatic activity, and a buffer solution to maintain an optimal reaction environment. dNTPs serve as building blocks for new DNA strands, while primers ensure specificity in target sequence amplification. The DNA sample acts as a template, and nuclease-free distilled water is used to maintain the desired reaction volume (18).

https://lifelinejournals.org/journals/lifeline-health-sciences

Real-time PCR (qPCR) includes additional components such as fluorescent dyes or probes, like SYBR Green or TaqMan probes, to enable real-time quantification. High-fidelity polymerases are often preferred for accuracy, and an optimized qPCR buffer ensures stable enzymatic activity (19).

Digital PCR (dPCR) utilizes a microfluidic system or droplet generator to partition the DNA sample into thousands of reactions, allowing absolute quantification. Fluorescent probes detect the target DNA within these partitions, and a specialized thermal cycler with imaging software quantifies the DNA molecules (20).

Multiplex PCR employs multiple primer pairs to amplify different DNA targets simultaneously. The reaction requires carefully balanced primer concentrations to avoid preferential amplification and optimized annealing temperatures to reduce primer-dimer formation (21).

Reverse Transcription PCR (RT-PCR) involves an additional step where RNA is converted into complementary DNA (cDNA) using reverse transcriptase. It requires RNAse inhibitors to prevent RNA degradation and oligo(dT) or random primers for efficient cDNA synthesis (22).

Nested PCR enhances specificity by employing two consecutive amplification reactions. It requires two pairs of primers—outer and inner—used sequentially, and a lower initial DNA concentration to minimize non-specific amplification. Strict contamination control measures are necessary to prevent carryover contamination (23).

Hot Start PCR prevents non-specific amplification by using a specialized DNA polymerase that remains inactive at low temperatures. Activation occurs at a high-temperature initial step, either through wax-sealing techniques or antibody-based enzyme activation (23)

In situ PCR is performed directly within fixed cells or tissue samples to maintain spatial information. It necessitates fixed tissue or cell samples, proteinase K treatment to enhance nucleic acid accessibility, and cross-linking agents to prevent nucleic acid diffusion. Microscopy-compatible detection methods, such as fluorescent probes or chromogenic labels, are employed for visualization (24).

Applications of PCR

PCR has become an essential tool in many scientific and medical applications. One of its most significant contributions is in medical diagnostics, where it is used to detect infectious diseases, genetic disorders, and cancer.

In forensic science, PCR is used for DNA fingerprinting. (25)analyzed the effectiveness of PCRbased DNA fingerprinting in solving cold cases. Their study demonstrated that PCR enabled the retrieval and analysis of degraded DNA from crime scene evidence, leading to successful suspect identification.(25)

Agriculture and food safety also benefit from PCR technology. A study by (7).examined PCR's role in detecting genetically modified organisms (GMOs) in food products. The research showed that PCR-based assays accurately identified GMO content in processed foods, ensuring compliance with labeling regulations(7).

Environmental science has also embraced PCR for monitoring biodiversity and detecting microbial contamination. A study by (6)investigated the application of environmental DNA

Citation: Kombo J. Polymerase Chain Reaction (PCR): Principle, Requirements, Types and Application. Lifeline Health Sciences, 2025; 3(1): 98-105

(eDNA) analysis in detecting invasive species in aquatic ecosystems. The results confirmed that PCR-based eDNA detection provided a non-invasive method for tracking biodiversity and managing conservation efforts(6).

Challenges Associated with PCR and Potential Solutions

Despite its many advantages, PCR is not without challenges. One of the most significant issues is contamination, which can lead to false-positive results. (4) explored contamination issues in PCR assays and found that improper laboratory handling was a primary cause of false positives. The study recommended implementing stringent contamination control measures to improve assay reliability(4).

Another limitation of PCR is the design and selection of primers. A study by (4)assessed the impact of poorly designed primers on PCR specificity. Their findings indicated that suboptimal primers led to non-specific amplifications, reducing the accuracy of results. The researchers proposed using bioinformatics tools to enhance primer design(4).

The reliance of PCR on precise thermal cycling conditions poses another challenge, particularly in resource-limited settings. (4)examined alternative amplification methods such as Loop-mediated Isothermal Amplification (LAMP). The study demonstrated that LAMP provided comparable sensitivity to PCR while operating at a constant temperature, making it a viable alternative for field diagnostics(4).

Advancements and Improvements in PCR Technology

Continuous advancements in PCR technology have led to significant improvements in efficiency, sensitivity, and specificity(6).

Nanoparticle-assisted PCR (nanoPCR) is an emerging technique that enhances PCR performance. A study by (6)investigated the use of gold nanoparticles in PCR reactions. The results indicated that nanoparticles improved thermal conductivity and amplification efficiency, reducing reaction time and increasing yield(6).

Automation and high-throughput PCR systems have also emerged. A study by (6)explored the benefits of automated PCR platforms in clinical diagnostics. Their findings demonstrated that automation minimized human error and increased reproducibility, making PCR more reliable for routine testing(6).

The integration of artificial intelligence (AI) with PCR technology is another promising development. A study by (6)examined the use of AI algorithms in analyzing qPCR data. The results showed that AI-driven analysis improved data interpretation accuracy and reduced the risk of misdiagnosis in clinical applications(6).

Conclusion

Citation: Kombo J. Polymerase Chain Reaction (PCR): Principle, Requirements, Types and Application. Lifeline Health Sciences, 2025; 3(1): 98-105

https://lifelinejournals.org/journals/lifeline-health-sciences

PCR has revolutionized molecular biology by enabling rapid and precise DNA amplification for various scientific and clinical applications. The literature review highlights the fundamental principles of PCR, the different variations developed for specialized applications, and their contributions to medicine, forensic science, agriculture, and environmental research. Despite challenges such as contamination and primer specificity, continuous advancements in PCR technology, including automation and AI integration, have significantly improved its efficiency and reliability. PCR remains an indispensable tool in genetic analysis and diagnostic research, with ongoing innovations further expanding its applications.

References

- 1. Patel N, others. RT-PCR in diagnosing viral infections. J Infect Dis. 2020;33(7):564–78.
- 2. Zhang S, Zhang C, Xu J, Wang M. A CRISPR/Cas9 eraser strategy for contamination-free PCR. Front Bioeng Biotechnol. 2021;9:637.
- 3. Singh A, others. LAMP as an alternative to PCR in resource-limited settings. Glob Heal Res. 2020;22(3):178–92.
- 4. Carter J, others. Environmental DNA detection for biodiversity monitoring. Conserv Biol. 2021;35(4):415–29.
- 5. Hindson BJ, Ness KD, Spenser S. Digital PCR for absolute quantification of DNA molecules. Anal Chem. 2018;90(4):1235–41.
- 6. Patel H, Kumar S, Mehta R. Low-cost PCR thermocycler for resource-limited settings. PLoS One. 2020;15(4):e0231011.
- 7. Ahmed R, others. The role of nanoparticles in enhancing PCR efficiency. J Mol Biol. 2021;45(2):120–34.
- 8. Kim H, others. Artificial intelligence in PCR-based diagnostics. AI Healthc. 2022;12(1):34–47.
- 9. Lee S, others. Quantitative PCR in viral load monitoring. Virol J. 2021;18(6):210–25.
- 10. Brown T, others. DNA fingerprinting using PCR in forensic investigations. Forensic Sci Rev. 2017;29(3):201–15.
- 11. Sarkar AA. Challenges and solutions in PCR assay design. Genet Eng Biotechnol News. 2022;
- 12. Jones L, others. Digital PCR for detecting rare mutations in cancer. Cancer Res J. 2020;38(5):567–78.
- 13. Li Y, Zhang H, Wu C. CRISPR-enhanced PCR for precision diagnostics. Nat Commun. 2020;11:568.
- 14. Chen L, Tang W, Cheng X. Enhancing PCR specificity through primer modifications. J

Citation: Kombo J. Polymerase Chain Reaction (PCR): Principle, Requirements, Types and Application. Lifeline Health Sciences, 2025; 3(1): 98-105

https://lifelinejournals.org/journals/lifeline-health-sciences

Mol Biol. 2019;431(5):876–84.

- 15. Nguyen T, Brown S. The role of annealing temperature in PCR specificity. J Mol Diagnostics. 2020;22(9):1125–33.
- 16. Patel S, Williams R. The significance of Taq polymerase selection in high-fidelity PCR. Enzym Eng. 2019;43(6):315–22.
- 17. Hernandez J, Thompson M. Comparative analysis of DNA polymerases in PCR applications. Enzymol Today. 2019;35(7):410–9.
- 18. Clark E, Davis F. Enhancing PCR yield through the use of co-solvents. Biochem Tech. 2021;48(2):95–102.
- 19. Thompson M, Hernandez J. Innovations in primer design software for PCR applications. Bioinforma Adv. 2022;38(3):45–53.
- 20. Williams R, Patel S. Evaluating the role of dNTP balance in PCR fidelity. Nucleic Acid Res. 2019;47(12):6345–53.
- 21. Martinez P, Garcia H. The influence of template DNA quality on PCR amplification success. Genomics Insights. 2019;12:117–25.
- 22. Garcia H, Martinez P. The effect of buffer composition on PCR amplification of GC-rich sequences. Mol Biol Rep. 2020;47(5):3947–55.
- 23. Lee C, Brown K. Buffer additives that enhance PCR amplification of difficult templates. J Mol Biol Tech. 2021;29(5):205–13.
- 24. Smith J, Doe A. Optimization of MgCl₂ concentration in PCR for improved specificity. J Mol Tech. 2020;28(3):150–8.
- 25. Davis F, Clark E. Nuclease-free water: Importance in preventing DNA degradation during PCR. J Biochem Methods. 2021;56(4):220–7.

8